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
or for BOWA handle with adjustable adapter bush and fingerswitch

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Medical Device Directive (Europe)

This device complies with the requirements of Directive 93/42 EEC concerning medical devices

1. Use and Operating

LAP electrodes are used for monopolar cutting and coagulation in combination with proper connecting cables and HF units. These products are more or less subject to wear depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. For this reason, please carry out a thorough visual examination and a function test prior to each application of these products. Cracks in the ceramic insulation may result from improper handling. In such cases, the product cannot be used any longer. A product must be replaced if it shows signs of wear or fails to function as described in these instructions. In this case, please contact the manufacturer or one of its authorized representatives.

Connection to HF unit in monopolar mode

In a first step, fix the neutral electrode ("return plate") to the patient, possibly on the upper arm or the thigh. Ensure that the skin is free from hair and grease. Apply conductive gel evenly over the neutral electrode. The current paths inside the patient's body should be short and must run diagonally. Current paths must never run diagonally across the body or across the thorax.

The patient must be insulated against all electrically conductive parts. Ground the operating table. Place the patient on a dry, electrically insulating layer. Prevent skin areas from touching each other and place dry gauze in between. Never place HF cables directly on the patient's skin. Never lay HF cables in loops. Connect instrument or handle with coagulating electrode and switch on the HF unit. Then connect the neutral electrode, the footswitch and the LAP electrode with the corresponding cable to the HF unit.

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In order to avoid interference with camera systems, do not lead HF cables in line with camera cables. Access to the surgical site is made by a trocar with a diameter of 5 mm. The maximum permissible high voltage of this electrode is 8500 V_{p-p}.

Always refer to the instruction manual of the HF unit and to its general information about electro-surgical procedures.

2. LAP electrodes with suction device

The models 710-xxx dispose of a suction tube and a universal hose connector for inner diameters of 5/6/7mm.

Sealing cap: 700-099

3. Cleaning and Disinfection

Cleaning*	
Manual / mechanic	Yes**
Ultrasound	Yes
Washing-up liquid	Yes
Machine	Yes
Disinfection*	
Manual	Yes
Machine	Yes
Sterilisation	
Autoclaving (134°C, 20 min)	Yes
Ethylene oxide	Yes
Formaldehyde	Yes
Plasma	No
Dry heat	No

* see separate recommendation sheet

** cleaning of suction tube with cleaning brush only

In the case of machine wash, the instructions of the machine producer must be respected.

4. Warnings

The products should only be used in combination with BOWA accessories and by clinically-trained professionals.

Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulating capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. Inspect for the correct position of the neutral electrode (in case of monopolar mode), correct attachment of all HF cables and plugs, correct activation of the current by the finger or footswitches, undamaged insulation of the HF cables and of the electrodes, and check if the distal end of the electrode is clean.

Before operating check if the signal transmission of the monitors works without interference (e.g. without noise); this check is done by pressing the yellow button (CUT) and the blue one (COAG).

Endogenous burns are those caused by high current density within the patient's tissue. Possible reasons may be that the patient is inadvertently positioned in contact with electrically conductive parts, or direct contact between skin areas and HF cables which may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by the heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons are among others: ignition of skin cleaning agents, ignition of narcotic gases etc.

Cardiac pacemakers may be damaged by electro-surgical current. Consult a cardiologist before the intervention. Never use electro-surgical current on outpatients fitted with pacemakers.

5. Disposal

The disposal of medical equipment, packaging material and accessories must be in accordance with relevant country-specific regulations and laws.