


ECOSAC

Semi-Automatic Range

- Cost effective tissue retrieval system
- Easy deployment
- Mouth or bottom first abdominal entry simplifies sac placement
- Large open mouth to facilitate tissue capture
- Secure drawstring closure

Product Code	Trocar (mm)	Mouth Diameter (mm)	Length (mm)	Volume (ml)	No. per box
EMP 70 ECO	10	57	100	100	10
EMP 100 ECO	10	70	144	200	10
EMP 140 ECO	12	102	175	550	10
EMP 170 ECO	12	170	170	1200	5
EMP 180 ECO	12	170	240	2000	5
New EMP 200 ECO	12	153	240	2500	5
New EMP 230 ECO	12	153	340	3100	5

*MDA is the sole distributor in Hong Kong, please contact us for more details.



MASTER E-SACS

Fully Automatic Tissue Retrieval System

- ① Automatically opening Sac
- ① Stronger arms reinforce mouth opening
- ① Drawstring ensures full closure
- ① Integral tail aids removal of Sac
- ① Accepts tissue sizes up to large Spleen, Ovary and Kidney

Product Code	Trocar (mm)	Mouth Diameter (mm)	Length (mm)	Volume (ml)	No. per box
EMP 80 WG	10	64	105	150	6
EMP 100 WG	10	64	150	200	6
EMP 120 WG	10	85	210	550	6
EMP 150 WG	12	110	170	950	3
EMP 160 WG	12	110	260	1500	3
EMP 170 WG	15	115	180	1000	3
EMP 180 WG	15	115	260	1800	3

Sac Size	Operation
70	Node biopsy, appendix, ectopic pregnancy, small gallbladder
80	Standard gallbladder, ovary, small uterus, fibroid, prostate, adrenal
100	Enlarged gallbladder, ovary, fibroid, adrenal, prostate, kidney
120	Ovarian cyst/tumour, large fibroid, uterus, kidney, spleen, bowel resection
140, 150 & 160	Malignant kidney, enlarged spleen, ovarian tumour, colon and lung resection
170 & 180	Large tumours of ovary, kidney, spleen and colon; lung and liver resection