

inside spread

210mm

210mm

208mm

A RANGE OF PRODUCTS TO MEET THE NEEDS OF MODERN SURGICAL PRACTICE

The introduction of the Chex™ CPH 34 and CPH 32 range of Circular Staplers for Stapled Hemorrhoidopexy has been received with the highest praise from Surgeons around the world. Undoubtedly...

"summa cum laude"

CPH 34 ACTUAL SIZE

21.60cm³ INTERNAL
HOUSING VOLUME


TOP VIEW



CPH 32 ACTUAL SIZE

17.17cm³ INTERNAL
HOUSING VOLUME

SIDE VIEW




ONLY PRODUCT THAT HAS 32 STAPLES
IN 2 ROWS TO ASSIST MORE EFFECTIVE
HEMOSTASIS

HIGHEST QUALITY SURGICAL STEEL
CUTTING BLADE

UNIQUE "CUTTING" MECHANISM WITH
CIRCUMFERENTIAL PUSH ASSEMBLY

RADICAL RE-DESIGN OF STAPLE HOUSING
GIVING SURGEON ACCURATE CONTROL
OF TISSUE DRAWN INTO HOUSING

LOWERED TISSUE STOP TO FACILITATE
SMOOTH CONTROLLABLE ENTRY OF
TISSUE INTO SUTURE HOUSING



A RANGE OF PRODUCTS THAT, QUITE SIMPLY, ARE NOW CONSIDERED THE BEST

When you put health first, you test everything as strenuously as possible. You test it yourself and have others test it for you, you ask questions, you go back if necessary, you incorporate what you've learned into the manufacturing process.

You devote all your energies to product quality and end-user cost reduction. Only then can you be satisfied that the products you produce are truly the best possible.

At Frankenman, we believe we're not so much in the business of making instruments but in the business of making surgeries more successful.

That's why we listen to doctors who've put our instruments to use. We read their reports. We invite them to think with us, to brainstorm, to help us develop new, innovative products.

That way we know we will be able to deliver the highest quality to the medical profession both today and way into the future.

CHEX™
HEALTH FIRST



ACCESSORIES



STAPLE INFORMATION

MODEL	CARTRIDGE COLOUR	STAPLE ROWS	OUTSIDE DIAMETER	CUTTING DIAMETER	STAPLE LEG LENGTH	STAPLE CROWN LENGTH	STAPLE COUNT	PER BOX
CPH 34	GREEN	2	34.5MM	26MM	4.2MM	3.8MM	32	3
CPH 32	WHITE	2	32.5MM	24MM	4.2MM	3.8MM	32	3

KEY FEATURES

- LARGE INTERNAL HOUSING CAPACITY TO ACCOMMODATE MORE TISSUE
- ADDITIONAL STAPLES TO PROVIDE EFFECTIVE HEMOSTASIS
- WING NUT CLOSING HANDLE FACILITATES EASE OF CLOSURE

outside spread

208mm

210mm

210mm

A RANGE OF PRODUCTS THAT FACILITATE MORE EFFECTIVE SURGERY



A commitment to the highest quality and an ability to exceed international standards is what makes Frankenman and our Chex range of products stand out. We put health first. This manifests itself in our total dedication to a zero defect outcome. Our manufacturing, production, research and development, and customer service facilities are the envy of the international medical device sector.

We source and use only the highest quality materials and components from around the world. This means that not only do we adhere to internationally recognized measures of quality, such as CE and ISO 13485 certification, we regularly exceed them. Our commitment to the ongoing training and development

of our staff is second to none and our design capabilities are of the highest standards throughout the industry culminating in the award of full ISO 9001: 2000 certification for our Quality Management systems.



PROFESSOR ULRICH BAUMGARTNER OF THE KREISKLINIK OTTOBUREN PRESENTS HIS EXPERIENCE OF THE USE OF THE CHEX CPH RANGE OF STAPLERS



PROFESSOR ULRICH
BAUMGARTNER

CHIEF OF THE DEPARTMENT OF
GENERAL, VISCERAL SURGERY AND
COLO-PROCTOLOGY SPECIALIST FOR
SURGERY AND VISCERAL SURGERY,
PROCTOLOGY (SPECIALIST AREAS OF
INTEREST: COLORECTAL SURGERY)

Professor Baumgartner is a recognised and experienced
trainer in Stapled Hemorrhoidopexy, STARR and other
colorectal procedures and has produced an excellent video
guide on the techniques he employs.

WHEN YOU PUT PATIENTS' HEALTH FIRST, REMARKABLE THINGS FOLLOW



State-of-the-art manufacturing facilities.

OUR RANGE OF PRODUCTS INCLUDE

SINGLE USE
RELOADABLE LINEAR
CUTTER STAPLER

SINGLE USE
RELOADABLE LINEAR
CUTTER STAPLER

SINGLE USE CURVED
INTRALUMINAL
CIRCULAR STAPLER



CHEX™
HEALTH FIRST

INFORMED DESIGN.
SCRUPULOUSLY APPLIED.
METICULOUSLY PRODUCED.

FOR MORE INFORMATION OR TO VIEW OUR ENTIRE PRODUCT RANGE PLEASE VISIT OUR WEBSITE.



SUITE B 13/F WING TAT COMMERCIAL BUILDING
121-125 WING LOK STREET SHEUNG WAN HONG KONG
T +852 3106 3035 F +852 3585 0519
E info@frankenman.hk

www.frankenman.hk
Chex and Chex Healthcare are registered trademarks of Frankenman International Limited



INFORMED DESIGN.
SCRUPULOUSLY APPLIED.
METICULOUSLY PRODUCED.



CIRCULAR STAPLERS FOR RECTAL PROLAPSE AND HEMORRHOIDS

CHEX™
HEALTH FIRST