




Kii®

Access Systems

First Entry System



Kii

Balloon Blunt Tip System

INNOVATIVE BLUNT TIP TECHNOLOGY



Suture Ties

Maximizes instrument reach and facilitates port site closure

Non-Latex Balloon

Provides stable anchoring throughout lengthy procedures while maximizing working space and visibility



Now available in a 5mm – the smallest hasson available!

Gel Cone

Enhances instrument articulation


Optical Blunt Tip Obturator

Facilitates direct visualization during trocar placement

Kii

Advanced Fixation Sleeve

SUPERIOR RETENTION



Retention Disk

Reduces migration of the trocar sleeve into the abdominal cavity




Non-Latex Balloon

Provides unmatched retention

Maximizes working space and visibility

Provides pressure on incision site, potentially reducing port-site bleeding

Maintains stability throughout procedures



Kii

Low Profile System

SLIM, LIGHTWEIGHT DESIGN

Low Profile Seal Housing

Optimizes external working space, minimizing interference with adjacent ports




Low Profile Kii Seal



Standard 5mm Kii Seal

Z-thread Sleeve with Extended Threads

Maintains retention at shorter depths, maximizing internal working space and visibility



Kii

Shielded Bladed System

SMOOTH, CONTROLLED ENTRY

Intuitive Switch Design

Provides clear indication of shield position



Deactivated



Activated

Kii Shielded Bladed System

Reorder No.		Size
<i>Advanced Fixation</i> <i>Z-thread</i>		
CFB03	CTB03	5 x 100mm
N/A	C0Q61	8 x 100mm
CFB33	CTB33	11 x 100mm
CFB73	CTB73	12 x 100mm

Kii

Optical Access System


UNPARALLELED OPTICAL CLARITY



Tissue Reapproximation



Initial
Withdrawal



Midway



Complete
Withdrawal

Kii Optical Access System

Reorder No.		Size
<i>Advanced Fixation</i>	<i>Z-thread</i>	
CFR03	CTR03	5 x 100mm
CFR01	CTR01	5 x 150mm
N/A	C0Q19	8 x 100mm
CFR33	CTR33	11 x 100mm
CFR31	CTR31	11 x 150mm
CFR73	CTR73	12 x 100mm
CFR71	CTR71	12 x 150mm
N/A	C0R37	15 x 100mm
N/A	C0R39	15 x 150mm



Kii Fios First Entry System

Reorder No.		Size
<i>Advanced Fixation</i>	<i>Z-thread</i>	
CFF05	N/A	5 x 75mm
CFF03	CTF03	5 x 100mm
CFF01	CTF01	5 x 150mm
CFF33	CTF33	11 x 100mm
CFF31	CTF31	11 x 150mm
CFF73	CTF73	12 x 100mm
CFF71	CTF71	12 x 150mm



Kii Advanced Fixation Sleeves

Reorder No.		Size
CFS03		5 x 75mm
CFS02		5 x 100mm
CFS01		5 x 150mm
CFS12		11 x 100mm
CFS11		11 x 150mm
CFS22		12 x 100mm



Kii Low Profile System

Reorder No.			Size
<i>Optical</i>	<i>Bladed</i>	<i>Sleeve</i>	
CTR21	CTB21	CTS28	5 x 55mm
CTR22	CTB22	CTS29	5 x 100mm
CTR23	CTB23	N/A	5 x 55mm, dual pack
CTR24	CTB24	N/A	5 x 100mm, dual pack



Kii Balloon Blunt Tip System

Reorder No.		Size
C0Q20		5 x 100mm
C0R47		12 x 100mm
C0R50		12 x 130mm

Visit www.appliedmedical.com/trocars
for more information