



Alexis[®]

Contained Extraction System



- Isolates specimen from surrounding cavity during tissue extraction
- Enables a minimally invasive approach
- Accommodates a wide variety of specimen sizes

Alexis[®] Contained Extraction System



Ring



Guard

Tether

Specimen
Containment Bag

Alexis Contained Extraction System

MODEL NO.	DESCRIPTION	BAG VOLUME	QUANTITY
GTK17	Alexis contained extraction system 17cm with Kii [®] balloon blunt tip system	6500 mL	1/KIT
GTM17	Alexis contained extraction system 17cm with GelPOINT [®] mini advanced access platform	6500 mL	1/KIT
GTG17	Alexis contained extraction system 17cm with GelPOINT advanced access platform	6500 mL	1/KIT
GTK14	Alexis contained extraction system 14cm with Kii balloon blunt tip system	3400 mL	1/KIT
GTM14	Alexis contained extraction system 14cm with GelPOINT mini advanced access platform	3400 mL	1/KIT
GTG14	Alexis contained extraction system 14cm with GelPOINT advanced access platform	3400 mL	1/KIT



<http://contact.appliedmedical.com/>

For more information
www.appliedmedical.com/alexiscs